

VD-8000A Biphasic Defibrillator



Technical specifications

■ Defibrillator

- Waveform: Biphasic
- Defibrillator type: Manual External Asynchronous Defibrillation.
- Charge Time: Less than 10 seconds to 360 Joules
- Charge Prompt Type: Voice and visual prompts.
- Energy Selection: Selectable at 0, 3, 5, 7, 10, 20, 30, 50, 100, 200, 300, 360 joules
- Standard Adult/Pediatric Paddles: Reusable external adult paddles (pediatric paddles integrated)

■ Battery Description

- Standard Type: Sealed lead battery FP1223C 12Vdc 2300mAh X2
- Standard Capacity: Up to 40 full-energy discharges.
- Low Battery Indicator: Bat Indicator on the panel of machine be on

■ Power

- AC-DC dual-use: Use the batteries, AC100V~240V, 50/60 Hz.
- Vehicle Voltage(Optional): DC 12V

■ Main unit:

- Portable
- Maximum power 600W
- Dimensions: L330xW220xH410
- Weight: 7.5kg