

# LVOT-005


## Electrical hydraulic Pet Operating Table



This is electrical hydraulic operating table use for bull and horse operation.  
It's made of stainless steel, equipped with IV pole, waste liquid tank and stainless steel tray.

### Main Spectification

Model	LVOT-005
Length	3700mm
Width	1200mm
Height	600-870mm
Table til	0-30°
Load Capacity	1500kg



Four detachable leg supports.  
Size: L700mm \* W400mm



Foot control pedal controls the lift and down of table top.



The hydraulic control system which can support 1500 kg weight.

Equipped with universal brake wheel.